

## Product datasheet for **TL303760**

### Girdin (CCDC88A) Human shRNA Plasmid Kit (Locus ID 55704)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Girdin (CCDC88A) Human shRNA Plasmid Kit (Locus ID 55704)
Locus ID:	55704
Synonyms:	APE; GIRDIN; GIV; GRDN; HkRP1; KIAA1212; PEHO; PEHOL
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CCDC88A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55704). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001135597</a> , <a href="#">NM_001254943</a> , <a href="#">NM_018084</a> , <a href="#">NM_018084.1</a> , <a href="#">NM_018084.2</a> , <a href="#">NM_018084.3</a> , <a href="#">NM_018084.4</a> , <a href="#">NM_001135597.1</a> , <a href="#">NM_001254943.1</a> , <a href="#">BC014905</a> , <a href="#">BC015539</a> , <a href="#">BC032683</a> , <a href="#">BC035848</a> , <a href="#">BC041160</a> , <a href="#">BC045829</a> , <a href="#">BC056885</a> , <a href="#">BC060037</a> , <a href="#">BC078168</a> , <a href="#">BC092515</a> , <a href="#">BC132736</a> , <a href="#">BC142646</a> , <a href="#">BC142700</a> , <a href="#">BC144320</a> , <a href="#">NM_001365480</a> , <a href="#">NM_018084.5</a> , <a href="#">NM_001135597.2</a> , <a href="#">NM_001254943.2</a>
UniProt ID:	<a href="#">Q3V6T2</a>
Summary:	This gene encodes a member of the Girdin family of coiled-coil domain containing proteins. The encoded protein is an actin-binding protein that is activated by the serine/threonine kinase Akt and plays a role in cytoskeleton remodeling and cell migration. The encoded protein also enhances Akt signaling by mediating phosphoinositide 3-kinase (PI3K)-dependent activation of Akt by growth factor receptor tyrosine kinases and G protein-coupled receptors. Increased expression of this gene and phosphorylation of the encoded protein may play a role in cancer metastasis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).