

## Product datasheet for **TL303797**

### **KCTD1 Human shRNA Plasmid Kit (Locus ID 284252)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	KCTD1 Human shRNA Plasmid Kit (Locus ID 284252)
<b>Locus ID:</b>	284252
<b>Synonyms:</b>	C18orf5
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	KCTD1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 284252). 5µg purified plasmid DNA per construct Non-effective 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<u><a href="#">NM_001136205</a></u> , <u><a href="#">NM_001142730</a></u> , <u><a href="#">NM_001258221</a></u> , <u><a href="#">NM_001258222</a></u> , <u><a href="#">NM_198991</a></u> , <u><a href="#">NM_001351443</a></u> , <u><a href="#">BC063652</a></u> , <u><a href="#">BC042371</a></u> , <u><a href="#">NM_001142730.3</a></u> , <u><a href="#">NM_001258222.3</a></u>
<b>Summary:</b>	This gene encodes a protein containing a BTB (Broad-complex, tramtrack and bric a brac), also known as a POZ (POxvirus and zinc finger) protein-protein interaction domain. The encoded protein negatively regulates the AP-2 family of transcription factors and the Wnt signaling pathway. A mechanism for the modulation of Wnt signaling has been proposed in which the encoded protein enhances ubiquitination and degradation of the beta-catenin protein. Mutations in this gene have been identified in Scalp-ear-nipple (SEN) syndrome. [provided by RefSeq, May 2017]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).