

Product datasheet for **TL303844**

KAZALD1 Human shRNA Plasmid Kit (Locus ID 81621)

Product data:

Product Type:	shRNA Plasmids
Product Name:	KAZALD1 Human shRNA Plasmid Kit (Locus ID 81621)
Locus ID:	81621
Synonyms:	BONO1; FKSG28; FKSG40; IGFBP-rP10
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	KAZALD1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 81621). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001319303 , NM_030929 , NR_135067 , NR_135068 , NM_030929.1 , NM_030929.2 , NM_030929.3 , NM_030929.4 , BC007758 , BC007758.2 , BC053611 , NM_030929.5
UniProt ID:	Q96I82
Summary:	This gene encodes a secreted member of the insulin growth factor-binding protein (IGFBP) superfamily. The protein contains an insulin growth factor-binding domain in its N-terminal region, a Kazal-type serine protease inhibitor and follistatin-like domain in its central region, and an immunoglobulin-like domain in its C-terminal region. Studies of the mouse ortholog suggest that this protein may function in bone development and bone regeneration. This gene is hypomethylated and over-expressed in high-grade glioma compared to low-grade glioma, and thus the hypomethylated gene may be associated with cell proliferation and the shorter survival of patients with high-grade glioma. It is also one of numerous genes found to be deleted in a novel 5.54 Mb interstitial deletion, which is associated with multiple congenital anomalies. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).