

Product datasheet for **TL303858V**

Junctional Adhesion Molecule 2 (JAM2) Human shRNA Lentiviral Particle (Locus ID 58494)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Junctional Adhesion Molecule 2 (JAM2) Human shRNA Lentiviral Particle (Locus ID 58494)
Locus ID:	58494
Synonyms:	C21orf43; CD322; IBGC8; JAM-B; JAMB; PRO245; VE-JAM; VEJAM
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	JAM2 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001270407 , NM_001270408 , NM_021219 , NR_072999 , NM_021219.1 , NM_021219.2 , NM_021219.3 , NM_001270407.1 , BC017779 , BC017779.1 , BM804795 , BM921497 , NM_021219.4 , NM_001270408.2
UniProt ID:	P57087
Summary:	This gene belongs to the immunoglobulin superfamily, and the junctional adhesion molecule (JAM) family. The protein encoded by this gene is a type I membrane protein that is localized at the tight junctions of both epithelial and endothelial cells. It acts as an adhesive ligand for interacting with a variety of immune cell types, and may play a role in lymphocyte homing to secondary lymphoid organs. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).