

## Product datasheet for **TL303910**

### INTS3 Human shRNA Plasmid Kit (Locus ID 65123)

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | shRNA Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Product Name:             | INTS3 Human shRNA Plasmid Kit (Locus ID 65123)                                                                                                                                                                                                                                                                                                                                                                                                 |
| Locus ID:                 | 65123                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Synonyms:                 | C1orf60; C1orf193; INT3; SOSS-A; SOSSA                                                                                                                                                                                                                                                                                                                                                                                                         |
| Vector:                   | pGFP-C-shLenti (TR30023)                                                                                                                                                                                                                                                                                                                                                                                                                       |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Format:                   | Lentiviral plasmids                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Components:               | INTS3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 65123).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.                                                                                                                                                                                                                   |
| RefSeq:                   | <a href="#">NM_001324475</a> , <a href="#">NM_023015</a> , <a href="#">NM_023015.1</a> , <a href="#">NM_023015.2</a> , <a href="#">NM_023015.3</a> , <a href="#">NM_023015.4</a> ,<br><a href="#">BC116458</a> , <a href="#">BC116458.1</a> , <a href="#">BC025254</a> , <a href="#">BC054513</a> , <a href="#">BC073985</a> , <a href="#">BC098431</a> , <a href="#">BC105092</a> , <a href="#">BC105094</a> ,<br><a href="#">NM_023015.5</a> |
| UniProt ID:               | <a href="#">Q68E01</a>                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Summary:                  | The protein encoded by this gene can form a complex with human single-strand DNA binding proteins 1 or 2 (hSSB1 and hSSB2) and other proteins to mediate genome stability and the DNA damage response. The encoded protein is also part of a multiprotein complex that interacts with the C-terminal domain of RNA polymerase II large subunit to help regulate processing of U1 and U2 small nuclear RNAs. [provided by RefSeq, May 2016]     |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .                                                                                             |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).