

## Product datasheet for **TL303968V**

### Cell adhesion molecule 2 (CADM2) Human shRNA Lentiviral Particle (Locus ID 253559)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cell adhesion molecule 2 (CADM2) Human shRNA Lentiviral Particle (Locus ID 253559)
Locus ID:	253559
Synonyms:	IGSF4D; Necl-3; NECL3; SynCAM 2; synCAM2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CADM2 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001167674</a> , <a href="#">NM_001167675</a> , <a href="#">NM_001256502</a> , <a href="#">NM_001256503</a> , <a href="#">NM_001256504</a> , <a href="#">NM_001256505</a> , <a href="#">NM_153184</a> , <a href="#">NM_153184.1</a> , <a href="#">NM_153184.2</a> , <a href="#">NM_153184.3</a> , <a href="#">NM_001167675.1</a> , <a href="#">NM_001167674.1</a> , <a href="#">NM_001256504.1</a> , <a href="#">NM_001256505.1</a> , <a href="#">NM_001256502.1</a> , <a href="#">NM_001256503.1</a> , <a href="#">BC105999</a> , <a href="#">BC040985</a> , <a href="#">BC063484</a> , <a href="#">BC112136</a> , <a href="#">BC112138</a> , <a href="#">BM674450</a> , <a href="#">NM_001167675.2</a> , <a href="#">NM_001167674.2</a> , <a href="#">NM_001256502.2</a> , <a href="#">NM_001256503.2</a> , <a href="#">NM_001256505.2</a> , <a href="#">NM_001256504.2</a>
UniProt ID:	<a href="#">Q8N3J6</a>
Summary:	This gene encodes a member of the synaptic cell adhesion molecule 1 (SynCAM) family which belongs to the immunoglobulin (Ig) superfamily. The encoded protein has three Ig-like domains and a cytosolic protein 4.1 binding site near the C-terminus. Proteins belonging to the protein 4.1 family crosslink spectrin and interact with other cytoskeletal proteins. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).