

Product datasheet for **TL304027**

HSH2D Human shRNA Plasmid Kit (Locus ID 84941)

Product data:

Product Type:	shRNA Plasmids
Product Name:	HSH2D Human shRNA Plasmid Kit (Locus ID 84941)
Locus ID:	84941
Synonyms:	ALX; HSH2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	HSH2D - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84941). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC025237 , NM_001291274 , NM_032855 , NR_111903 , NR_111904 , NM_001352265 , NM_001352266 , NM_032855.1 , NM_032855.2 , NM_032855.3 , NM_001291274.1 , BC025237.2 , BC016826 , NM_001369808 , NR_163150 , NR_163152 , NR_163153 , NR_163155 , NM_001369809 , NR_163151 , NR_163154 , NR_163156
UniProt ID:	Q96JZ2
Summary:	T-cell activation requires 2 signals: recognition of antigen by the T-cell receptor (see TCR; MIM 186880) and a costimulatory signal provided primarily by CD28 (MIM 186760) in naive T cells. HSH2 is a target of both of these signaling pathways (Greene et al., 2003 [PubMed 12960172]).[supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).