

Product datasheet for **TL304096**

HIVEP3 Human shRNA Plasmid Kit (Locus ID 59269)

Product data:

Product Type:	shRNA Plasmids
Product Name:	HIVEP3 Human shRNA Plasmid Kit (Locus ID 59269)
Locus ID:	59269
Synonyms:	KBP-1; KBP1; KRC; Schnurri-3; SHN3; ZAS3; ZNF40C
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	HIVEP3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 59269). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001127714 , NM_024503 , NR_038260 , NR_038261 , NM_024503.1 , NM_024503.2 , NM_024503.3 , NM_024503.4 , NM_001127714.1 , NM_001127714.2 , BC015242 , BC040942 , BC046120 , BC062792 , BC080552 , BC152563 , BC172754 , NM_001127714.3 , NM_024503.5
UniProt ID:	Q5T1R4
Summary:	This gene encodes a member of the human immunodeficiency virus type 1 enhancer-binding protein family. Members of this protein family contain multiple zinc finger and acid-rich (ZAS) domains and serine-threonine rich regions. This protein acts as a transcription factor and is able to regulate nuclear factor kappaB-mediated transcription by binding the kappaB motif in target genes. This protein also binds the recombination signal sequence that flanks the V, D, and J regions of immunoglobulin and T-cell receptors. Alternate splicing results in both coding and non-coding transcript variants. [provided by RefSeq, Sep 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).