

Product datasheet for **TL304403V**

GAMT Human shRNA Lentiviral Particle (Locus ID 2593)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	GAMT Human shRNA Lentiviral Particle (Locus ID 2593)
Locus ID:	2593
Synonyms:	CCDS2; HEL-S-20; PIG2; TP53I2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	GAMT - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_000156 , NM_138924 , NM_000156.3 , NM_000156.4 , NM_000156.5 , NM_138924.1 , NM_138924.2 , BC016760 , BC016760.2 , BC017936 , BM994704 , NM_000156.6 , NM_138924.3
UniProt ID:	Q14353
Summary:	The protein encoded by this gene is a methyltransferase that converts guanidoacetate to creatine, using S-adenosylmethionine as the methyl donor. Defects in this gene have been implicated in neurologic syndromes and muscular hypotonia, probably due to creatine deficiency and accumulation of guanidinoacetate in the brain of affected individuals. Two transcript variants encoding different isoforms have been described for this gene. Pseudogenes of this gene are found on chromosomes 2 and 13. [provided by RefSeq, Feb 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).