

Product datasheet for **TL304466**

FREM1 Human shRNA Plasmid Kit (Locus ID 158326)

Product data:

Product Type:	shRNA Plasmids
Product Name:	FREM1 Human shRNA Plasmid Kit (Locus ID 158326)
Locus ID:	158326
Synonyms:	BNAR; C9orf143; C9orf145; C9orf154; MOTA; TILRR; TRIGNO2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FREM1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 158326). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC031064 , NM_001177704 , NM_144966 , NM_144966.1 , NM_144966.3 , NM_144966.4 , NM_144966.5 , NM_001177704.1 , BC031064.1 , BC172226 , NM_001370060 , NM_001370061 , NM_001370065 , NR_163239 , NR_163240 , NR_163241 , NM_001370058 , NM_001370063 , NR_163238 , NR_163242 , NM_001177704.3 , NM_144966.7
UniProt ID:	Q5H8C1
Summary:	This gene encodes a basement membrane protein that may play a role in craniofacial and renal development. Mutations in this gene have been associated with bifid nose with or without anorectal and renal anomalies. Alternatively spliced transcript variants encoding different isoforms have been described. PubMed ID 19940113 describes one such variant that initiates transcription within a distinct, internal exon; the resulting shorter isoform (named Toll-like/interleukin-1 receptor regulator, TILRR) is suggested to be a co-receptor of the interleukin 1 receptor family and may regulate receptor function and Toll-like receptor/interleukin 1 receptor signal transduction, contributing to the control of inflammatory response activation. [provided by RefSeq, Apr 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).