

## Product datasheet for **TL304476**

### FOXI1 Human shRNA Plasmid Kit (Locus ID 2299)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	FOXI1 Human shRNA Plasmid Kit (Locus ID 2299)
Locus ID:	2299
Synonyms:	FKH10; FKHL10; FREAC-6; FREAC6; HFH-3; HFH3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FOXI1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 2299). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_012188</a> , <a href="#">NM_144769</a> , <a href="#">NM_144769.1</a> , <a href="#">NM_144769.2</a> , <a href="#">NM_012188.1</a> , <a href="#">NM_012188.2</a> , <a href="#">NM_012188.3</a> , <a href="#">NM_012188.4</a> , <a href="#">BC029778</a> , <a href="#">BC029778.1</a> , <a href="#">NM_012188.5</a>
UniProt ID:	<a href="#">Q12951</a>
Summary:	This gene belongs to the forkhead family of transcription factors, which is characterized by a distinct forkhead domain. This gene may play an important role in the development of the cochlea and vestibulum, as well as in embryogenesis. The encoded protein has been found to be required for the transcription of four subunits of a proton pump found in the inner ear, the kidney, and the epididymis. Mutations in this gene have been associated with deafness, autosomal recessive 4. [provided by RefSeq, Jan 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).