

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TL304489

Fibronectin (FN1) Human shRNA Plasmid Kit (Locus ID 2335)

Product data:

| Product Type: | shRNA Plasmids |
|------------------------------|---|
| Product Name: | Fibronectin (FN1) Human shRNA Plasmid Kit (Locus ID 2335) |
| Locus ID: | 2335 |
| Synonyms: | CIG; ED-B; FINC; FN; FNZ; GFND; GFND2; LETS; MSF; SMDCF |
| Vector: | pGFP-C-shLenti (TR30023) |
| E. coli Selection: | Chloramphenicol (34 ug/ml) |
| Mammalian Cell Selection: | Puromycin |
| Format: | Lentiviral plasmids |
| Components: | FN1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 2335). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free. |
| RefSeq: | NM 001306129, NM 001306130, NM 001306131, NM 001306132, NM 002026, NM 054034, NM 212474, NM 212475, NM 212476, NM 212478, NM 212482, NM 054034.1, NM 054034.2, NM 212482.1, NM 212482.2, NM 212474.1, NM 212474.2, NM 212476.1, NM 212476.2, NM 212478.1, NM 212478.2, NM 002026.1, NM 002026.2, NM 002026.3, NM 212475.1, BC117176, BC005858, BC016875, BC078656, BC100030, BC143754, BC143763, NM 001365517, NM 001365520, NM 001365518, NM 001365519, NM 001365521, NM 001365522, NM 001365523, NM 001365524, NM 212482.3, NM 054034.3, NM 212474.3 |
| UniProt ID: | <u>P02751</u> |
| Summary: | This gene encodes fibronectin, a glycoprotein present in a soluble dimeric form in plasma, and in a dimeric or multimeric form at the cell surface and in extracellular matrix. The encoded preproprotein is proteolytically processed to generate the mature protein. Fibronectin is involved in cell adhesion and migration processes including embryogenesis, wound healing, blood coagulation, host defense, and metastasis. The gene has three regions subject to alternative splicing, with the potential to produce 20 different transcript variants, at least one of which encodes an isoform that undergoes proteolytic processing. The full-length nature of some variants has not been determined. [provided by RefSeq, Jan 2016] |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

| | Fibronectin (FN1) Human shRNA Plasmid Kit (Locus ID 2335) – TL304489 |
|----------------------------|--|
| shRNA Design: | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u> . If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u> . |
| Performance Guaranteed: | OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples. |
| | For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred). |

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US