

Product datasheet for **TL304492**

FMN2 Human shRNA Plasmid Kit (Locus ID 56776)

Product data:

Product Type:	shRNA Plasmids
Product Name:	FMN2 Human shRNA Plasmid Kit (Locus ID 56776)
Locus ID:	56776
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FMN2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 56776). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001305424 , NM_020066 , NM_001348094 , NM_020066.4 , NM_020066.1 , NM_020066.3 , BC014364 , BC112335 , BC112361
UniProt ID:	Q9NZ56
Summary:	This gene is a member of the formin homology protein family. The encoded protein is thought to have essential roles in organization of the actin cytoskeleton and in cell polarity. This protein mediates the formation of an actin mesh that positions the spindle during oogenesis and also regulates the formation of actin filaments in the nucleus. This protein also forms a perinuclear actin/focal-adhesion system that regulates the shape and position of the nucleus during cell migration. Mutations in this gene have been associated with infertility and also with an autosomal recessive form of intellectual disability (MRT47). Alternatively spliced transcript variants have been identified. [provided by RefSeq, Jul 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).