

## Product datasheet for **TL304560**

### **FBXO18 Human shRNA Plasmid Kit (Locus ID 84893)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	FBXO18 Human shRNA Plasmid Kit (Locus ID 84893)
<b>Locus ID:</b>	84893
<b>Synonyms:</b>	Fbx18; FBXO18; hFBH1
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	FBXO18 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84893). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<u><a href="#">NM_001258452</a></u> , <u><a href="#">NM_001258453</a></u> , <u><a href="#">NM_032807</a></u> , <u><a href="#">NM_178150</a></u> , <u><a href="#">NM_032807.1</a></u> , <u><a href="#">NM_032807.2</a></u> , <u><a href="#">NM_032807.3</a></u> , <u><a href="#">NM_032807.4</a></u> , <u><a href="#">NM_178150.1</a></u> , <u><a href="#">NM_178150.2</a></u> , <u><a href="#">NM_001258453.1</a></u> , <u><a href="#">NM_001258452.1</a></u> , <u><a href="#">BC006430</a></u> , <u><a href="#">BC012762</a></u> , <u><a href="#">BC020266</a></u> , <u><a href="#">BC032674</a></u> , <u><a href="#">BC080629</a></u> , <u><a href="#">BC110884</a></u> , <u><a href="#">BC113375</a></u> , <u><a href="#">BC113377</a></u> , <u><a href="#">NM_178150.3</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q8NFZ0</a></u>
<b>Summary:</b>	This gene encodes a member of the F-box protein family, members of which are characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into three classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbx class. It contains an F-box motif and seven conserved helicase motifs, and has both DNA-dependent ATPase and DNA unwinding activities. Alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).