

Product datasheet for **TL304993**

DIABLO Human shRNA Plasmid Kit (Locus ID 56616)

Product data:

Product Type:	shRNA Plasmids
Product Name:	DIABLO Human shRNA Plasmid Kit (Locus ID 56616)
Locus ID:	56616
Synonyms:	DFNA64; SMAC
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	DIABLO - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 56616). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001278302 , NM_001278303 , NM_001278304 , NM_001278342 , NM_019887 , NM_138929 , NM_138930 , NR_024600 , NR_024601 , NM_019887.1 , NM_019887.2 , NM_019887.3 , NM_019887.4 , NM_019887.5 , NM_138930.1 , NM_138930.2 , NM_138930.3 , NM_138929.1 , NM_138929.2 , NM_138929.3 , NM_001278302.1 , NM_001278304.1 , NM_001278342.1 , BC011909 , BC011909.2 , BC004417 , BC046209 , BC095468 , BM746689 , NM_019887.6 , NM_001278304.2
UniProt ID:	Q9NR28
Summary:	This gene encodes an inhibitor of apoptosis protein (IAP)-binding protein. The encoded mitochondrial protein enters the cytosol when cells undergo apoptosis, and allows activation of caspases by binding to inhibitor of apoptosis proteins. Overexpression of the encoded protein sensitizes tumor cells to apoptosis. A mutation in this gene is associated with young-adult onset of nonsyndromic deafness-64. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).