

Product datasheet for **TL305038**

DENND1C Human shRNA Plasmid Kit (Locus ID 79958)

Product data:

Product Type:	shRNA Plasmids
Product Name:	DENND1C Human shRNA Plasmid Kit (Locus ID 79958)
Locus ID:	79958
Synonyms:	FAM31C
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	DENND1C - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 79958). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001290331 , NM_024898 , NM_024898.1 , NM_024898.2 , NM_024898.3 , NM_001290331.1 , BC033437 , BC033437.1 , NM_001290331.2 , NM_024898.4
UniProt ID:	Q8IV53
Summary:	The protein encoded by this gene functions as a guanine nucleotide exchange factor for the early endosomal small GTPase RAB35, which regulates endosomal membrane trafficking and is involved in actin polymerization. The encoded protein activates RAB35 by promoting the exchange of RAB35-bound GDP for GTP. This gene may play a role in linking RAB35 activation with the clathrin machinery. [provided by RefSeq, May 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).