

Product datasheet for **TL305201V**

CSAD Human shRNA Lentiviral Particle (Locus ID 51380)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	CSAD Human shRNA Lentiviral Particle (Locus ID 51380)
Locus ID:	51380
Synonyms:	CSD; PCAP
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CSAD - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001244705 , NM_001244706 , NM_015989 , NM_015989.1 , NM_015989.2 , NM_015989.3 , NM_015989.4 , NM_001244706.1 , NM_001244705.1 , BC018042 , BC047466 , BC052249 , BC098278 , BC098342 , BC099717 , BC105918 , NM_015989.5 , NM_001244705.2
UniProt ID:	Q9Y600
Summary:	This gene encodes a member of the group 2 decarboxylase family. A similar protein in rodents plays a role in multiple biological processes as the rate-limiting enzyme in taurine biosynthesis, catalyzing the decarboxylation of cysteinesulfinate to hypotaurine. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).