

## Product datasheet for **TL305228V**

### CRBN Human shRNA Lentiviral Particle (Locus ID 51185)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	CRBN Human shRNA Lentiviral Particle (Locus ID 51185)
Locus ID:	51185
Synonyms:	MRT2; MRT2A
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CRBN - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001173482</a> , <a href="#">NM_016302</a> , <a href="#">NM_016302.1</a> , <a href="#">NM_016302.2</a> , <a href="#">NM_016302.3</a> , <a href="#">NM_001173482.1</a> , <a href="#">BC067811</a> , <a href="#">BC067811.1</a> , <a href="#">BC017419</a> , <a href="#">BC059171</a> , <a href="#">NM_016302.4</a>
UniProt ID:	<a href="#">Q96SW2</a>
Summary:	This gene encodes a protein related to the Lon protease protein family. In rodents and other mammals this gene product is found in the cytoplasm localized with a calcium channel membrane protein, and is thought to play a role in brain development. Mutations in this gene are associated with autosomal recessive nonsyndromic cognitive disability. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**

GFP signal was observed under microscope at 48 hours after transduction of TL305228B virus into HEK293 cells. TL305228B virus was prepared using lenti-shRNA TL305228B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL305228C] virus into HEK293 cells. [TL305228C] virus was prepared using lenti-shRNA [TL305228C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL305228D] virus into HEK293 cells. [TL305228D] virus was prepared using lenti-shRNA [TL305228D] and [TR30037] packaging kit.