

## Product datasheet for **TL305319**

### CNOT7 Human shRNA Plasmid Kit (Locus ID 29883)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CNOT7 Human shRNA Plasmid Kit (Locus ID 29883)
Locus ID:	29883
Synonyms:	CAF-1; CAF1; Caf1a; hCAF-1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CNOT7 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 29883). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_013354</a> , <a href="#">NM_054026</a> , <a href="#">NM_001322087</a> , <a href="#">NM_001322088</a> , <a href="#">NM_001322089</a> , <a href="#">NM_001322090</a> , <a href="#">NM_001322091</a> , <a href="#">NM_001322092</a> , <a href="#">NM_001322093</a> , <a href="#">NM_001322094</a> , <a href="#">NM_001322095</a> , <a href="#">NM_001322096</a> , <a href="#">NM_001322097</a> , <a href="#">NM_001322098</a> , <a href="#">NM_001322099</a> , <a href="#">NM_001322100</a> , <a href="#">NM_013354.1</a> , <a href="#">NM_013354.2</a> , <a href="#">NM_013354.3</a> , <a href="#">NM_013354.4</a> , <a href="#">NM_013354.5</a> , <a href="#">NM_013354.6</a> , <a href="#">NM_054026.1</a> , <a href="#">NM_054026.2</a> , <a href="#">NM_054026.3</a> , <a href="#">BC060852</a> , <a href="#">BC060852.1</a> , <a href="#">BC007315</a> , <a href="#">BC070187</a> , <a href="#">BM468592</a> , <a href="#">NM_013354.7</a>
UniProt ID:	<a href="#">Q9UIV1</a>
Summary:	The protein encoded by this gene binds to an anti-proliferative protein, B-cell translocation protein 1, which negatively regulates cell proliferation. Binding of the two proteins, which is driven by phosphorylation of the anti-proliferative protein, causes signaling events in cell division that lead to changes in cell proliferation associated with cell-cell contact. The encoded protein downregulates the innate immune response and therefore provides a therapeutic target for enhancing its antimicrobial activity against foreign agents. Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 1 and X. [provided by RefSeq, Apr 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).