

Product datasheet for **TL305354**

Dectin 1 (CLEC7A) Human shRNA Plasmid Kit (Locus ID 64581)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Dectin 1 (CLEC7A) Human shRNA Plasmid Kit (Locus ID 64581)
Locus ID:	64581
Synonyms:	BGR; CANDF4; CD369; CLECSF12; DECTIN1; SCARE2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CLEC7A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 64581). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_022570 , NM_197947 , NM_197948 , NM_197949 , NM_197950 , NM_197951 , NM_197952 , NM_197953 , NM_197954 , NR_125336 , NM_197949.1 , NM_197949.2 , NM_197950.1 , NM_197950.2 , NM_022570.1 , NM_022570.2 , NM_022570.3 , NM_022570.4 , NM_197947.1 , NM_197947.2 , NM_197948.1 , NM_197948.2 , NM_197954.1 , NM_197954.2 , NM_197951.1 , NM_197952.1 , BC013385 , BC071746 , BC093829 , BC093831 , NM_197949.3 , NM_022570.5 , NM_197948.3 , NM_197947.3 , NM_197950.3 , NM_197954.3
UniProt ID:	Q9BXN2
Summary:	This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. The encoded glycoprotein is a small type II membrane receptor with an extracellular C-type lectin-like domain fold and a cytoplasmic domain with an immunoreceptor tyrosine-based activation motif. It functions as a pattern-recognition receptor that recognizes a variety of beta-1,3-linked and beta-1,6-linked glucans from fungi and plants, and in this way plays a role in innate immune response. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region. [provided by RefSeq, Jul 2008]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).