

Product datasheet for **TL305463V**

CDYL Human shRNA Lentiviral Particle (Locus ID 9425)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	CDYL Human shRNA Lentiviral Particle (Locus ID 9425)
Locus ID:	9425
Synonyms:	CDYL1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CDYL - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001143970 , NM_001143971 , NM_004824 , NM_170751 , NM_170752 , NR_026590 , NM_004824.1 , NM_004824.2 , NM_004824.3 , NM_170752.1 , NM_001143970.1 , NM_001143971.1 , NM_170751.1 , BC108725 , BC043397 , BC061516 , BC119682 , NM_001368126 , NM_001368125 , NM_001368127 , NM_004824.4
UniProt ID:	Q9Y232
Summary:	Chromodomain Y is a primate-specific Y-chromosomal gene family expressed exclusively in the testis and implicated in infertility. Although the Y-linked genes are testis-specific, this autosomal gene is ubiquitously expressed. The Y-linked genes arose by retrotransposition of an mRNA from this gene, followed by amplification of the retroposed gene. Proteins encoded by this gene superfamily possess a chromodomain, a motif implicated in chromatin binding and gene suppression, and a catalytic domain believed to be involved in histone acetylation. Multiple proteins are encoded by transcript variants of this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).