

## Product datasheet for **TL305475**

### CDIPT Human shRNA Plasmid Kit (Locus ID 10423)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CDIPT Human shRNA Plasmid Kit (Locus ID 10423)
Locus ID:	10423
Synonyms:	PIS; PIS1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CDIPT - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10423). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001286585</a> , <a href="#">NM_001286586</a> , <a href="#">NM_006319</a> , <a href="#">NM_145752</a> , <a href="#">NM_006319.1</a> , <a href="#">NM_006319.2</a> , <a href="#">NM_006319.3</a> , <a href="#">NM_006319.4</a> , <a href="#">NM_001286586.1</a> , <a href="#">NM_001286585.1</a> , <a href="#">NM_145752.1</a> , <a href="#">BC001444</a> , <a href="#">BC001444.2</a> , <a href="#">NM_001286586.2</a> , <a href="#">NM_006319.5</a> , <a href="#">NM_001286585.2</a>
UniProt ID:	<a href="#">O14735</a>
Summary:	Phosphatidylinositol breakdown products are ubiquitous second messengers that function downstream of many G protein-coupled receptors and tyrosine kinases regulating cell growth, calcium metabolism, and protein kinase C activity. Two enzymes, CDP-diacylglycerol synthase and phosphatidylinositol synthase, are involved in the biosynthesis of phosphatidylinositol. Phosphatidylinositol synthase, a member of the CDP-alcohol phosphatidyl transferase class-I family, is an integral membrane protein found on the cytoplasmic side of the endoplasmic reticulum and the Golgi apparatus. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).