

## Product datasheet for **TL305497V**

### CDC42 Human shRNA Lentiviral Particle (Locus ID 998)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	CDC42 Human shRNA Lentiviral Particle (Locus ID 998)
Locus ID:	998
Synonyms:	CDC42Hs; G25K; TKS
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CDC42 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001039802</a> , <a href="#">NM_001791</a> , <a href="#">NM_044472</a> , <a href="#">NM_001791.1</a> , <a href="#">NM_001791.2</a> , <a href="#">NM_001791.3</a> , <a href="#">NM_044472.1</a> , <a href="#">NM_044472.2</a> , <a href="#">NM_001039802.1</a> , <a href="#">BC018266</a> , <a href="#">BC018266.1</a> , <a href="#">BC002711</a> , <a href="#">BC003682</a> , <a href="#">BC018622</a> , <a href="#">BM684474</a> , <a href="#">NM_001791.4</a> , <a href="#">NM_044472.3</a> , <a href="#">NM_001039802.2</a>
UniProt ID:	<a href="#">P60953</a>
Summary:	The protein encoded by this gene is a small GTPase of the Rho-subfamily, which regulates signaling pathways that control diverse cellular functions including cell morphology, migration, endocytosis and cell cycle progression. This protein is highly similar to <i>Saccharomyces cerevisiae</i> Cdc 42, and is able to complement the yeast <i>cdc42-1</i> mutant. The product of oncogene <i>Dbl</i> was reported to specifically catalyze the dissociation of GDP from this protein. This protein could regulate actin polymerization through its direct binding to Neural Wiskott-Aldrich syndrome protein (N-WASP), which subsequently activates Arp2/3 complex. Alternative splicing of this gene results in multiple transcript variants. Pseudogenes of this gene have been identified on chromosomes 3, 4, 5, 7, 8 and 20. [provided by RefSeq, Apr 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL305497A virus into HEK293 cells. TL305497A virus was prepared using lenti-shRNA TL305497A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL305497B virus into HEK293 cells. TL305497B virus was prepared using lenti-shRNA TL305497B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL305497C] virus into HEK293 cells. [TL305497C] virus was prepared using lenti-shRNA [TL305497C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL305497D] virus into HEK293 cells. [TL305497D] virus was prepared using lenti-shRNA [TL305497D] and [TR30037] packaging kit.