

## Product datasheet for TL305527

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CCR11 (ACKR4) Human shRNA Plasmid Kit (Locus ID 51554)

**Product data:** 

**Product Type:** shRNA Plasmids

**Product Name:** CCR11 (ACKR4) Human shRNA Plasmid Kit (Locus ID 51554)

Locus ID:

CC-CKR-11; CCBP2; CCR-11; CCR10; CCR11; CCRL1; CCX-CKR; CCX CKR; CKR-11; PPR1; VSHK1 Synonyms:

pGFP-C-shLenti (TR30023) Vector: E. coli Selection: Chloramphenicol (34 ug/ml)

Mammalian Cell

Selection:

Puromycin

Format:

Lentiviral plasmids

Components: ACKR4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51554).

5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.

NM 016557, NM 178445, NM 016557.1, NM 016557.2, NM 016557.3, NM 178445.1, RefSeq:

NM 178445.2, BC069438, BC069438.1, BC095501, NM 016557.4

**UniProt ID:** O9NPB9

**Summary:** The protein encoded by this gene is a member of the G protein-coupled receptor family, and

> is a receptor for C-C type chemokines. This receptor has been shown to bind dendritic celland T cell-activated chemokines including CCL19/ELC, CCL21/SLC, and CCL25/TECK. A

pseudogene of this gene is found on chromosome 6. Alternatively spliced transcript variants

encoding the same protein have been described. [provided by RefSeq, Jul 2013]

These shRNA constructs were designed against multiple splice variants at this gene locus. To shRNA Design:

> be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.







## Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).