

## Product datasheet for **TL305536V**

### Cyclin D1 (CCND1) Human shRNA Lentiviral Particle (Locus ID 595)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cyclin D1 (CCND1) Human shRNA Lentiviral Particle (Locus ID 595)
Locus ID:	595
Synonyms:	BCL1; D11S287E; PRAD1; U21B31
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CCND1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_053056</a> , <a href="#">NM_053056.1</a> , <a href="#">NM_053056.2</a> , <a href="#">BC023620</a> , <a href="#">BC023620.2</a> , <a href="#">BC000076</a> , <a href="#">BC001501</a> , <a href="#">BC014078</a> , <a href="#">BC025302</a> , <a href="#">BM796500</a>
UniProt ID:	<a href="#">P24385</a>
Summary:	The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance throughout the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin forms a complex with and functions as a regulatory subunit of CDK4 or CDK6, whose activity is required for cell cycle G1/S transition. This protein has been shown to interact with tumor suppressor protein Rb and the expression of this gene is regulated positively by Rb. Mutations, amplification and overexpression of this gene, which alters cell cycle progression, are observed frequently in a variety of human cancers. [provided by RefSeq, Dec 2019]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).