

## Product datasheet for **TL305694V**

### CA12 Human shRNA Lentiviral Particle (Locus ID 771)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	CA12 Human shRNA Lentiviral Particle (Locus ID 771)
Locus ID:	771
Synonyms:	CA-XII; CAXII; HsT18816; T18816
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CA12 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001218</a> , <a href="#">NM_001293642</a> , <a href="#">NM_206925</a> , <a href="#">NR_135511</a> , <a href="#">NM_206925.1</a> , <a href="#">NM_206925.2</a> , <a href="#">NM_001218.1</a> , <a href="#">NM_001218.2</a> , <a href="#">NM_001218.3</a> , <a href="#">NM_001218.4</a> , <a href="#">NM_001293642.1</a> , <a href="#">BC011691</a> , <a href="#">BC011691.2</a> , <a href="#">BC023981</a> , <a href="#">BC023981.1</a> , <a href="#">BC000278</a> , <a href="#">BC001012</a> , <a href="#">BC087838</a> , <a href="#">BM994285</a> , <a href="#">NM_001218.5</a>
UniProt ID:	<a href="#">O43570</a>
Summary:	Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. This gene product is a type I membrane protein that is highly expressed in normal tissues, such as kidney, colon and pancreas, and has been found to be overexpressed in 10% of clear cell renal carcinomas. Three transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jun 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).