

Product datasheet for **TL305881**

MEMO1 Human shRNA Plasmid Kit (Locus ID 51072)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MEMO1 Human shRNA Plasmid Kit (Locus ID 51072)
Locus ID:	51072
Synonyms:	C2orf4; CGI-27; MEMO; NS5ATP7
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MEMO1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51072). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001137602 , NM_001301833 , NM_001301852 , NM_015955 , NR_126032 , NR_126034 , NM_015955.1 , NM_015955.2 , NM_015955.3 , NM_001137602.1 , NM_001137602.2 , NM_001301852.1 , NM_001301833.1 , BC018733 , BC036262 , BC070046 , BC094681 , NM_015955.4 , NM_001301852.3 , NM_001301833.3 , NM_001137602.3
UniProt ID:	Q9Y316
Summary:	May control cell migration by relaying extracellular chemotactic signals to the microtubule cytoskeleton. Mediator of ERBB2 signaling. The MEMO1-RHOA-DIAPH1 signaling pathway plays an important role in ERBB2-dependent stabilization of microtubules at the cell cortex. It controls the localization of APC and CLASP2 to the cell membrane, via the regulation of GSK3B activity. In turn, membrane-bound APC allows the localization of the MACF1 to the cell membrane, which is required for microtubule capture and stabilization. Is required for breast carcinoma cell migration.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).