

Product datasheet for **TL306307**

PAMM (FAM213A) Human shRNA Plasmid Kit (Locus ID 84293)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PAMM (FAM213A) Human shRNA Plasmid Kit (Locus ID 84293)
Locus ID:	84293
Synonyms:	Adrx; C10orf58; FAM213A; PAMM
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FAM213A - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84293). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001136047 , NM_001243778 , NM_001243779 , NM_001243780 , NM_001243781 , NM_001243782 , NM_032333 , NR_024572 , NM_032333.1 , NM_032333.2 , NM_032333.3 , NM_032333.4 , NM_001243781.1 , NM_001243782.1 , NM_001243778.1 , NM_001243779.1 , NM_001243780.1 , BC005871 , BC024308 , BC045777 , BC071938 , NM_001243781.2 , NM_001243782.2 , NM_032333.5
UniProt ID:	Q9BRX8
Summary:	Involved in redox regulation of the cell (PubMed:26438880, PubMed:19951071). Acts as an antioxidant (PubMed:19951071, PubMed:26438880). Inhibits TNFSF11-induced NFKB1 and JUN activation and osteoclast differentiation (PubMed:19951071). May affect bone resorption and help to maintain bone mass (PubMed:19951071). Acts as a negative regulator of macrophage-mediated inflammation by inhibiting macrophage production of inflammatory cytokines, probably through suppression of the MAPK signaling pathway (PubMed:26438880). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).