

Product datasheet for **TL306375V**

BRF1 Human shRNA Lentiviral Particle (Locus ID 2972)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	BRF1 Human shRNA Lentiviral Particle (Locus ID 2972)
Locus ID:	2972
Synonyms:	BRF; BRF-1; CFDS; GTF3B; hBRF; HEL-S-76p; TAF3B2; TAF3C; TAFIII90; TF3B90; TFIIIB90
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	BRF1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001242786 , NM_001242787 , NM_001242788 , NM_001242789 , NM_001242790 , NM_001519 , NM_145685 , NM_145696 , NM_145685.2 , NM_001519.1 , NM_001519.2 , NM_001519.3 , NM_001242790.1 , NM_001242789.1 , NM_001242787.1 , NM_001242786.1 , NM_001242788.1 , NM_145696.1 , BC086856 , BC016743 , BC071637 , BC103859 , BC103860 , BC103861 , BC133655 , NM_001519.4 , NM_001242790.2 , NM_001242788.2
UniProt ID:	Q92994
Summary:	This gene encodes one of the three subunits of the RNA polymerase III transcription factor complex. This complex plays a central role in transcription initiation by RNA polymerase III on genes encoding tRNA, 5S rRNA, and other small structural RNAs. The gene product belongs to the TF2B family. Several alternatively spliced variants encoding different isoforms, that function at different promoters transcribed by RNA polymerase III, have been identified. [provided by RefSeq, Jun 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).