

## Product datasheet for **TL306423**

### **BCL2L10 Human shRNA Plasmid Kit (Locus ID 10017)**

#### **Product data:**

|                                  |  |
|----------------------------------|--|
| <b>Product Type:</b>             | shRNA Plasmids   |
| <b>Product Name:</b>             | BCL2L10 Human shRNA Plasmid Kit (Locus ID 10017)   |
| <b>Locus ID:</b>                 | 10017  |
| <b>Synonyms:</b>                 | BCL-B; bcl2-L-10; Boo; Diva  |
| <b>Vector:</b>                   | pGFP-C-shLenti (TR30023)   |
| <b>E. coli Selection:</b>        | Chloramphenicol (34 ug/ml)   |
| <b>Mammalian Cell Selection:</b> | Puromycin  |
| <b>Format:</b>                   | Lentiviral plasmids  |
| <b>Components:</b>               | BCL2L10 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10017). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.  |
| <b>RefSeq:</b>                   | <a href="#">NM_001306168</a> , <a href="#">NM_020396</a> , <a href="#">NM_020396.1</a> , <a href="#">NM_020396.2</a> , <a href="#">NM_020396.3</a> , <a href="#">BC104442</a> , <a href="#">BC093826</a> , <a href="#">BC093828</a> , <a href="#">BC104443</a> , <a href="#">BC143227</a>  |
| <b>UniProt ID:</b>               | <a href="#">Q9HD36</a>   |
| <b>Summary:</b>                  | The protein encoded by this gene belongs to the BCL-2 protein family. BCL-2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. The protein encoded by this gene contains conserved BH4, BH1 and BH2 domains. This protein can interact with other members of BCL-2 protein family including BCL2, BCL2L1/BCL-X(L), and BAX. Overexpression of this gene has been shown to suppress cell apoptosis possibly through the prevention of cytochrome C release from the mitochondria, and thus activating caspase-3 activation. The mouse counterpart of this protein is found to interact with Apaf1 and forms a protein complex with Caspase 9, which suggests the involvement of this protein in APAF1 and CASPASE 9 related apoptotic pathway. [provided by RefSeq, Jul 2008] |
| <b>shRNA Design:</b>             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).