

## Product datasheet for **TL306442V**

### **BAT5 (ABHD16A) Human shRNA Lentiviral Particle (Locus ID 7920)**

#### **Product data:**

Product Type:	shRNA Lentiviral Particles
Product Name:	BAT5 (ABHD16A) Human shRNA Lentiviral Particle (Locus ID 7920)
Locus ID:	7920
Synonyms:	BAT5; D6S82E; hBAT5; NG26; PP199
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	BAT5 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001177515</a> , <a href="#">NM_021160</a> , <a href="#">NR_033488</a> , <a href="#">NR_033489</a> , <a href="#">NM_021160.1</a> , <a href="#">NM_021160.2</a> , <a href="#">NM_001177515.1</a> , <a href="#">BC031839</a> , <a href="#">NM_021160.3</a> , <a href="#">NM_001177515.2</a>
UniProt ID:	<a href="#">O95870</a>
Summary:	A cluster of genes, BAT1-BAT5, has been localized in the vicinity of the genes for tumor necrosis factor alpha and tumor necrosis factor beta. These genes are all within the human major histocompatibility complex class III region. The protein encoded by this gene is thought to be involved in some aspects of immunity. Alternatively spliced transcript variants have been described. [provided by RefSeq, Apr 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).