

Product datasheet for **TL306494**

ATP8A2 Human shRNA Plasmid Kit (Locus ID 51761)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ATP8A2 Human shRNA Plasmid Kit (Locus ID 51761)
Locus ID:	51761
Synonyms:	ATP; ATP1B; CAMRQ4; IB; ML-1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ATP8A2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51761). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_016529 , NM_001313741 , NM_016529.1 , NM_016529.2 , NM_016529.3 , NM_016529.4 , NM_016529.5 , BC144228 , BC144229 , BC156471 , BC172535 , NM_016529.6
UniProt ID:	Q9NTI2
Summary:	The protein encoded by this gene is a member of the P4 ATPase family of proteins, which are thought to be involved in a process called lipid flipping, whereby phospholipids are translocated inwards from the exoplasmic leaflet to the cytosolic leaflet of the cell membrane, which aids in generating and maintaining asymmetry in membrane lipids. This protein is predicted to contain an E1 E2 ATPase, a haloacid dehalogenase-like hydrolase (HAD) domain, and multiple transmembrane domains. Associations between this protein and cell cycle control protein 50A are important for translocation of phosphatidylserine across membranes. Mutations in this gene have been associated with a syndrome (CAMRQ4) characterized by cerebellar ataxia and cognitive disabilities. In addition, a translocation breakpoint within this gene was observed in an individual with neurological dysfunction. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).