

## Product datasheet for **TL306532V**

### ASPM Human shRNA Lentiviral Particle (Locus ID 259266)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	ASPM Human shRNA Lentiviral Particle (Locus ID 259266)
Locus ID:	259266
Synonyms:	ASP; Calmbp1; MCPH5
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	ASPM - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, $>10^7$ TU/ml.
RefSeq:	<a href="#">NM_001206846</a> , <a href="#">NM_018136</a> , <a href="#">NM_018136.1</a> , <a href="#">NM_018136.2</a> , <a href="#">NM_018136.3</a> , <a href="#">NM_018136.4</a> , <a href="#">NM_001206846.1</a> , <a href="#">BC015396</a> , <a href="#">BC034607</a> , <a href="#">BC040439</a> , <a href="#">NM_018136.5</a>
UniProt ID:	<a href="#">Q8IZT6</a>
Summary:	This gene is the human ortholog of the Drosophila melanogaster 'abnormal spindle' gene (asp), which is essential for normal mitotic spindle function in embryonic neuroblasts. Studies in mouse also suggest a role of this gene in mitotic spindle regulation, with a preferential role in regulating neurogenesis. Mutations in this gene are associated with microcephaly primary type 5. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).