

Product datasheet for **TL306670**

ANKS1B Human shRNA Plasmid Kit (Locus ID 56899)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ANKS1B Human shRNA Plasmid Kit (Locus ID 56899)
Locus ID:	56899
Synonyms:	AIDA; AIDA-1; ANKS2; cajalin-2; EB-1; EB1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ANKS1B - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 56899). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001204065 , NM_001204066 , NM_001204067 , NM_001204068 , NM_001204069 , NM_001204070 , NM_001204079 , NM_001204080 , NM_001204081 , NM_020140 , NM_152788 , NM_181670 , NM_001352185 , NM_001352186 , NM_001352187 , NM_001352188 , NM_001352189 , NM_001352190 , NM_001352191 , NM_001352192 , NM_001352193 , NM_001352194 , NM_001352195 , NM_001352196 , NM_001352197 , NM_001352198 , NM_001352199 , NM_001352200 , NM_001352201 , NM_001352202 , NM_001352203 , NM_001352204 , NM_001352205 , NM_001352206 , NM_001352207 , NM_001352208 , NM_001352209 , NM_001352210 , NM_001352211 , NM_001352212 , NM_001352213 , NM_001352214 , NM_001352216 , NM_001352217 , NM_001352218 , NM_001352219 , NM_001352220 , NM_001352221 , NM_001352222 , NM_001352223 , NM_001352224 , NM_001352225 , NM_152788.1 , NM_152788.2 , NM_152788.3 , NM_152788.4 , NM_181670.1 , NM_181670.2 , NM_181670.3 , NM_020140.1 , NM_020140.2 , NM_020140.3 , NM_001204081.1 , NM_001204080.1 , NM_001204065.1 , NM_001204079.1 , NM_001204066.1 , NM_001204067.1 , NM_001204070.1 , NM_001204069.1 , NM_001204068.1 , BC026313 , BC068451 , BC091512 , BC142669 , BC150204 , BC160005 , NM_001204070.2 , NM_001204066.2 , NM_181670.4 , NM_001204068.2 , NM_001204065.2 , NM_001204080.2 , NM_001204079.2 , NM_001204069.2
UniProt ID:	Q7Z6G8



[View online »](#)

- Summary:** This gene encodes a multi-domain protein that is predominantly expressed in brain and testis. This protein interacts with amyloid beta protein precursor (AbetaPP) and may have a role in normal brain development, and in the pathogenesis of Alzheimer's disease. Expression of this gene has been shown to be elevated in patients with pre-B cell acute lymphocytic leukemia associated with t(1;19) translocation. Alternatively spliced transcript variants encoding different isoforms (some with different subcellular localization, PMID:15004329) have been described for this gene. [provided by RefSeq, Aug 2011]
- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).