

## Product datasheet for **TL306709**

### ANKHD1 Human shRNA Plasmid Kit (Locus ID 54882)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ANKHD1 Human shRNA Plasmid Kit (Locus ID 54882)
Locus ID:	54882
Synonyms:	MASK; MASK1; PP2500; VBARP
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ANKHD1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 54882). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001197030</a> , <a href="#">NM_017747</a> , <a href="#">NM_017978</a> , <a href="#">NM_024668</a> , <a href="#">NM_020690.2</a> , <a href="#">NM_017978.1</a> , <a href="#">NM_017978.2</a> , <a href="#">NM_024668.1</a> , <a href="#">NM_024668.2</a> , <a href="#">NM_024668.3</a> , <a href="#">NM_017747.1</a> , <a href="#">NM_017747.2</a> , <a href="#">NM_001197030.1</a> , <a href="#">BC004457</a> , <a href="#">BC009420</a> , <a href="#">BC009909</a> , <a href="#">BC040231</a> , <a href="#">BC117677</a> , <a href="#">BC117678</a> , <a href="#">BC127127</a> , <a href="#">BC150486</a> , <a href="#">BC172415</a> , <a href="#">NM_017978.3</a> , <a href="#">NM_024668.4</a> , <a href="#">NM_017747.3</a>
UniProt ID:	<a href="#">Q8IWZ3</a>
Summary:	This gene encodes a protein with multiple ankyrin repeat domains and a single KH-domain. The protein is thought to function as a scaffolding protein, and it may be involved in the regulation of caspases and thereby play an antiapoptotic role in cell survival. Alternative splicing results in multiple transcript variants, one of which generates a fusion transcript (MASK-BP3) with the downstream eIF4E-binding protein 3 (EIF4EBP3) gene, resulting in a protein comprised of the ANKHD1 sequence for the majority of the protein and a different C-terminus due to an alternate reading frame for the EIF4EBP3 segments. [provided by RefSeq, Sep 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).