

Product datasheet for **TL306849**

ADAM15 Human shRNA Plasmid Kit (Locus ID 8751)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ADAM15 Human shRNA Plasmid Kit (Locus ID 8751)
Locus ID:	8751
Synonyms:	MDC15
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ADAM15 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 8751). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001261464 , NM_001261465 , NM_001261466 , NM_003815 , NM_207191 , NM_207194 , NM_207195 , NM_207196 , NM_207197 , NR_048577 , NR_048578 , NR_048579 , NM_003815.1 , NM_003815.2 , NM_003815.3 , NM_003815.4 , NM_207191.1 , NM_207191.2 , NM_207194.1 , NM_207194.2 , NM_207195.1 , NM_207195.2 , NM_207197.1 , NM_207197.2 , NM_207196.1 , NM_207196.2 , NM_001261466.1 , NM_001261465.1 , NM_001261464.1 , BC014566 , BC014566.2 , BM971880 , NM_001261465.2 , NM_001261466.2 , NM_207196.3 , NM_207197.3 , NM_001261464.2 , NM_207191.3 , NM_207195.3 , NM_207194.3 , NM_003815.5
UniProt ID:	Q13444
Summary:	The protein encoded by this gene is a member of the ADAM (a disintegrin and metalloproteinase) protein family. ADAM family members are type I transmembrane glycoproteins known to be involved in cell adhesion and proteolytic ectodomain processing of cytokines and adhesion molecules. This protein contains multiple functional domains including a zinc-binding metalloprotease domain, a disintegrin-like domain, as well as a EGF-like domain. Through its disintegrin-like domain, this protein specifically interacts with the integrin beta chain, beta 3. It also interacts with Src family protein-tyrosine kinases in a phosphorylation-dependent manner, suggesting that this protein may function in cell-cell adhesion as well as in cellular signaling. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).