

## Product datasheet for **TL307109**

### THNSL2 Human shRNA Plasmid Kit (Locus ID 55258)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	THNSL2 Human shRNA Plasmid Kit (Locus ID 55258)
Locus ID:	55258
Synonyms:	SOFAT; THS2; TSH2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	THNSL2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55258). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001244676</a> , <a href="#">NM_001244678</a> , <a href="#">NM_018271</a> , <a href="#">NM_018271.1</a> , <a href="#">NM_018271.2</a> , <a href="#">NM_018271.3</a> , <a href="#">NM_018271.4</a> , <a href="#">NM_001244678.1</a> , <a href="#">NM_001244676.1</a> , <a href="#">BC035315</a> , <a href="#">BC047758</a> , <a href="#">BC050699</a> , <a href="#">BC064423</a> , <a href="#">NM_018271.5</a>
UniProt ID:	<a href="#">Q86YI6</a>
Summary:	This gene encodes a threonine synthase-like protein. A similar enzyme in mouse can catalyze the degradation of O-phospho-homoserine to a-ketobutyrate, phosphate, and ammonia. This protein also has phospho-lyase activity on both gamma and beta phosphorylated substrates. In mouse an alternatively spliced form of this protein has been shown to act as a cytokine and can induce the production of the inflammatory cytokine IL6 in osteoblasts. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).