

## Product datasheet for **TL307137**

### Ubiquitin (UBB) Human shRNA Plasmid Kit (Locus ID 7314)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Ubiquitin (UBB) Human shRNA Plasmid Kit (Locus ID 7314)
Locus ID:	7314
Synonyms:	HEL-S-50
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	UBB - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 7314). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001281716</a> , <a href="#">NM_001281717</a> , <a href="#">NM_001281718</a> , <a href="#">NM_001281719</a> , <a href="#">NM_001281720</a> , <a href="#">NM_018955</a> , <a href="#">NM_018955.1</a> , <a href="#">NM_018955.2</a> , <a href="#">NM_018955.3</a> , <a href="#">NM_001281716.1</a> , <a href="#">NM_001281717.1</a> , <a href="#">NM_001281718.1</a> , <a href="#">NM_001281719.1</a> , <a href="#">NM_001281720.1</a> , <a href="#">BC000379</a> , <a href="#">BC000379.2</a> , <a href="#">BC009301</a> , <a href="#">BC015127</a> , <a href="#">BC026301</a> , <a href="#">BC031027</a> , <a href="#">BC038999</a> , <a href="#">BC046123</a> , <a href="#">BM554254</a> , <a href="#">BM662083</a> , <a href="#">NM_018955.4</a>
UniProt ID:	<a href="#">P0CG47</a>
Summary:	This gene encodes ubiquitin, one of the most conserved proteins known. Ubiquitin has a major role in targeting cellular proteins for degradation by the 26S proteasome. It is also involved in the maintenance of chromatin structure, the regulation of gene expression, and the stress response. Ubiquitin is synthesized as a precursor protein consisting of either polyubiquitin chains or a single ubiquitin moiety fused to an unrelated protein. This gene consists of three direct repeats of the ubiquitin coding sequence with no spacer sequence. Consequently, the protein is expressed as a polyubiquitin precursor with a final amino acid after the last repeat. An aberrant form of this protein has been detected in patients with Alzheimer's disease and Down syndrome. Pseudogenes of this gene are located on chromosomes 1, 2, 13, and 17. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]



[View online »](#)

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).