

Product datasheet for **TL307170**

LYRM1 Human shRNA Plasmid Kit (Locus ID 57149)

Product data:

Product Type:	shRNA Plasmids
Product Name:	LYRM1 Human shRNA Plasmid Kit (Locus ID 57149)
Locus ID:	57149
Synonyms:	A211C6.1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	LYRM1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 57149). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001128301 , NM_001128302 , NM_001302835 , NM_001302836 , NM_020424 , NR_126516 , NR_126520 , NR_126521 , NR_126527 , NR_126528 , NM_020424.1 , NM_020424.2 , NM_020424.4 , NM_001128301.1 , NM_001128301.2 , NM_001128302.1 , NM_001128302.2 , NM_001302835.1 , NM_001302836.1 , BC017039 , BC017039.1 , BM546065 , NM_001369632 , NM_001369633 , NM_001369635 , NM_001369636 , NM_001369638 , NM_001369639 , NM_001369641 , NM_001369642 , NM_001369644 , NR_161425 , NR_161426 , NR_161429 , NR_161430 , NM_001369631 , NM_001369634 , NM_001369637 , NM_001369645 , NR_161427 , NR_161431 , NR_161432 , NM_001302835.2 , NM_001302836.2 , NM_001128302.3 , NM_001128301.3 , NM_020424.5
UniProt ID:	O43325
Summary:	The protein encoded by this gene belongs to the mitochondrial leucine/tyrosine/arginine motif family of proteins. Proteins of this family are short polypeptides that contain a leucine/tyrosine/arginine motif near the N-terminus. This gene is widely expressed with high levels in omental adipose tissue of obese individuals. In adipose tissue, the protein is localized to the nucleus where it promotes preadipocyte proliferation and lowers the rate of apoptosis to regulate adipose tissue homeostasis. Overexpression of this gene in adipocytes causes abnormal mitochondrial morphology and mitochondrial dysfunction. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2014]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).