

Product datasheet for **TL307183V**

Raptor (RPTOR) Human shRNA Lentiviral Particle (Locus ID 57521)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Raptor (RPTOR) Human shRNA Lentiviral Particle (Locus ID 57521)
Locus ID:	57521
Synonyms:	KOG1; Mip1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	RPTOR - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001163034 , NM_020761 , NM_020761.1 , NM_020761.2 , NM_001163034.1 , BC014502 , BC025180 , BC033258 , BC033699 , BC041079 , BC064515 , BC073898 , BC136652 , BC136654 , NM_001163034.2 , NM_020761.3
UniProt ID:	Q8N122
Summary:	This gene encodes a component of a signaling pathway that regulates cell growth in response to nutrient and insulin levels. The encoded protein forms a stoichiometric complex with the mTOR kinase, and also associates with eukaryotic initiation factor 4E-binding protein-1 and ribosomal protein S6 kinase. The protein positively regulates the downstream effector ribosomal protein S6 kinase, and negatively regulates the mTOR kinase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:

GFP signal was observed under microscope at 48 hours after transduction of TL307183B virus into HEK293 cells. TL307183B virus was prepared using lenti-shRNA TL307183B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL307183D] virus into HEK293 cells. [TL307183D] virus was prepared using lenti-shRNA [TL307183D] and [TR30037] packaging kit.