

## Product datasheet for **TL307421**

### MICALL1 Human shRNA Plasmid Kit (Locus ID 85377)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	MICALL1 Human shRNA Plasmid Kit (Locus ID 85377)
Locus ID:	85377
Synonyms:	MICAL-L1; MIRAB13
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MICALL1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 85377). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_033386</a> , <a href="#">NM_033386.1</a> , <a href="#">NM_033386.2</a> , <a href="#">NM_033386.3</a> , <a href="#">BC001090</a> , <a href="#">BC082243</a> , <a href="#">BC142605</a> , <a href="#">BC143050</a> , <a href="#">BC146383</a> , <a href="#">BC148804</a> , <a href="#">NM_033386.4</a>
UniProt ID:	<a href="#">Q8N3F8</a>
Summary:	Probable lipid-binding protein with higher affinity for phosphatidic acid, a lipid enriched in recycling endosome membranes. On endosome membranes, may act as a downstream effector of Rab proteins recruiting cytosolic proteins to regulate membrane tubulation. May be involved in a late step of receptor-mediated endocytosis regulating for instance endocytosed-EGF receptor trafficking. Alternatively, may regulate slow endocytic recycling of endocytosed proteins back to the plasma membrane. May indirectly play a role in neurite outgrowth.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).