

Product datasheet for **TL307576**

SEPT6 Human shRNA Plasmid Kit (Locus ID 23157)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SEPT6 Human shRNA Plasmid Kit (Locus ID 23157)
Locus ID:	23157
Synonyms:	SEP2; SEPT2; SEPT6
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SEPT6 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 23157). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_015129 , NM_145799 , NM_145800 , NM_145802 , NM_145802.1 , NM_145802.2 , NM_145802.3 , NM_015129.1 , NM_015129.2 , NM_015129.3 , NM_015129.4 , NM_015129.5 , NM_145799.1 , NM_145799.2 , NM_145800.1 , NM_145800.2 , NM_145800.3 , BC011922 , BC036240 , BC069231 , BC069231.1 , BC009291 , NM_145800.4 , NM_145802.4 , NM_015129.6 , NM_145799.4
UniProt ID:	Q14141
Summary:	This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).