

## Product datasheet for **TL307858**

### ARS2 (SRRT) Human shRNA Plasmid Kit (Locus ID 51593)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ARS2 (SRRT) Human shRNA Plasmid Kit (Locus ID 51593)
Locus ID:	51593
Synonyms:	ARS2; ASR2; serrate
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SRRT - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51593). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001128852</a> , <a href="#">NM_001128853</a> , <a href="#">NM_001128854</a> , <a href="#">NM_015908</a> , <a href="#">NM_182800</a> , <a href="#">NM_015908.1</a> , <a href="#">NM_015908.3</a> , <a href="#">NM_015908.4</a> , <a href="#">NM_015908.5</a> , <a href="#">NM_001128854.1</a> , <a href="#">NM_001128853.1</a> , <a href="#">NM_001128852.1</a> , <a href="#">NM_182800.2</a> , <a href="#">BC109117</a> , <a href="#">BC000082</a> , <a href="#">BC069249</a> , <a href="#">BC109116</a> , <a href="#">NM_015908.6</a> , <a href="#">NM_001128852.2</a> , <a href="#">NM_001128854.2</a> , <a href="#">NM_001128853.2</a>
UniProt ID:	<a href="#">Q9BXP5</a>
Summary:	Acts as a mediator between the cap-binding complex (CBC) and the primary microRNAs (miRNAs) processing machinery during cell proliferation. Contributes to the stability and delivery of capped primary miRNA transcripts to the primary miRNA processing complex containing DGCR8 and DROSHA, thereby playing a role in RNA-mediated gene silencing (RNAi) by miRNAs. Binds capped RNAs (m7GpppG-capped RNA); however interaction is probably mediated via its interaction with NCBP1/CBP80 component of the CBC complex. Involved in cell cycle progression at S phase. Does not directly confer arsenite resistance but rather modulates arsenic sensitivity. Independently of its activity on miRNAs, necessary and sufficient to promote neural stem cell self-renewal. Does so by directly binding SOX2 promoter and positively regulating its transcription (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).