

Product datasheet for **TL307876**

STING (TMEM173) Human shRNA Plasmid Kit (Locus ID 340061)

Product data:

Product Type:	shRNA Plasmids
Product Name:	STING (TMEM173) Human shRNA Plasmid Kit (Locus ID 340061)
Locus ID:	340061
Synonyms:	ERIS; hMITA; hSTING; MITA; MPYS; NET23; SAVI; STING; STING-beta; TMEM173
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TMEM173 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 340061). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001301738 , NM_198282 , NM_198282.1 , NM_198282.2 , NM_198282.3 , NM_001301738.1 , BC047779 , BC047779.1 , BM149879 , BM981219 , NM_001367258 , NM_001301738.2 , NM_198282.4
UniProt ID:	Q86WV6
Summary:	This gene encodes a five transmembrane protein that functions as a major regulator of the innate immune response to viral and bacterial infections. The encoded protein is a pattern recognition receptor that detects cytosolic nucleic acids and transmits signals that activate type I interferon responses. The encoded protein has also been shown to play a role in apoptotic signaling by associating with type II major histocompatibility complex. Mutations in this gene are the cause of infantile-onset STING-associated vasculopathy. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL307876A virus into HEK293 cells. TL307876A virus was prepared using lenti-shRNA TL307876A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL307876B virus into HEK293 cells. TL307876B virus was prepared using lenti-shRNA TL307876B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL307876C] virus into HEK293 cells. [TL307876C] virus was prepared using lenti-shRNA [TL307876C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL307876D] virus into HEK293 cells. [TL307876D] virus was prepared using lenti-shRNA [TL307876D] and [TR30037] packaging kit.