

Product datasheet for **TL308084**

MPEG1 Human shRNA Plasmid Kit (Locus ID 219972)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MPEG1 Human shRNA Plasmid Kit (Locus ID 219972)
Locus ID:	219972
Synonyms:	Mpg-1; MPG1; MPS1; P-2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MPEG1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 219972). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001039396 , NM_001039396.1 , BC112230 , BC104997
UniProt ID:	Q2M385
Summary:	Plays a key role in the innate immune response following bacterial infection by inserting into the bacterial surface to form pores (By similarity). By breaching the surface of phagocytosed bacteria, allows antimicrobial effectors to enter the bacterial periplasmic space and degrade bacterial proteins such as superoxide dismutase sodC which contributes to bacterial virulence (By similarity). Shows antibacterial activity against a wide spectrum of Gram-positive, Gram-negative and acid-fast bacteria (PubMed:23753625, PubMed:26402460, PubMed:30609079). Reduces the viability of the intracytosolic pathogen L.monocytogenes by inhibiting acidification of the phagocytic vacuole of host cells which restricts bacterial translocation from the vacuole to the cytosol (By similarity). Required for the antibacterial activity of reactive oxygen species and nitric oxide (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).