

Product datasheet for **TL308278V**

ZFX Human shRNA Lentiviral Particle (Locus ID 7543)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	ZFX Human shRNA Lentiviral Particle (Locus ID 7543)
Locus ID:	7543
Synonyms:	ZNF926
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	ZFX - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001178084 , NM_001178085 , NM_001178086 , NM_001178095 , NM_003410 , NM_001330327 , NM_003410.1 , NM_003410.2 , NM_003410.3 , NM_001178095.1 , NM_001178086.1 , NM_001178084.1 , NM_001178085.1 , BC030756 , BC136308 , BC136309 , NM_003410.4 , NM_001178084.2 , NM_001178095.2
UniProt ID:	P17010
Summary:	This gene on the X chromosome is structurally similar to a related gene on the Y chromosome. It encodes a member of the krueppel C2H2-type zinc-finger protein family. The full-length protein contains an acidic transcriptional activation domain (AD), a nuclear localization sequence (NLS) and a DNA binding domain (DBD) consisting of 13 C2H2-type zinc fingers. Studies in mouse embryonic and adult hematopoietic stem cells showed that this gene was required as a transcriptional regulator for self-renewal of both stem cell types, but it was dispensable for growth and differentiation of their progeny. Multiple alternatively spliced transcript variants encoding different isoforms have been identified, but the full-length nature of some variants has not been determined. [provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).