

## Product datasheet for **TL308537**

### UBA52 Human shRNA Plasmid Kit (Locus ID 7311)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	UBA52 Human shRNA Plasmid Kit (Locus ID 7311)
Locus ID:	7311
Synonyms:	CEP52; HUBCEP52; L40; RPL40
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	UBA52 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 7311). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001033930</a> , <a href="#">NM_003333</a> , <a href="#">NM_001321017</a> , <a href="#">NM_001321018</a> , <a href="#">NM_001321019</a> , <a href="#">NM_001321020</a> , <a href="#">NM_001321021</a> , <a href="#">NM_001321022</a> , <a href="#">NM_003333.1</a> , <a href="#">NM_003333.2</a> , <a href="#">NM_003333.3</a> , <a href="#">NM_003333.4</a> , <a href="#">NM_001033930.1</a> , <a href="#">NM_001033930.2</a> , <a href="#">BC101830</a> , <a href="#">BC101830.1</a> , <a href="#">BC101832</a> , <a href="#">BC143668</a> , <a href="#">NM_001033930.3</a> , <a href="#">NM_003333.5</a>
UniProt ID:	<a href="#">P62987</a>
Summary:	Ubiquitin is a highly conserved nuclear and cytoplasmic protein that has a major role in targeting cellular proteins for degradation by the 26S proteasome. It is also involved in the maintenance of chromatin structure, the regulation of gene expression, and the stress response. Ubiquitin is synthesized as a precursor protein consisting of either polyubiquitin chains or a single ubiquitin moiety fused to an unrelated protein. This gene encodes a fusion protein consisting of ubiquitin at the N terminus and ribosomal protein L40 at the C terminus, a C-terminal extension protein (CEP). Multiple processed pseudogenes derived from this gene are present in the genome. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).