

Product datasheet for TL308538

OriGene Technologies, Inc.9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

UAP1 Human shRNA Plasmid Kit (Locus ID 6675)

Product data:

Product Type: shRNA Plasmids

Product Name: UAP1 Human shRNA Plasmid Kit (Locus ID 6675)

Locus ID: 6675

Synonyms: AGX; AGX1; AGX2; SPAG2

Vector: pGFP-C-shLenti (TR30023)

E. coli Selection: Chloramphenicol (34 ug/ml)

Mammalian Cell

Selection:

Puromycin

Format: Lentiviral plasmids

Components: UAP1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6675).

5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.

RefSeq: NM 003115, NM 001324113, NM 001324114, NM 001324115, NM 001324116,

NM 001324117, NM 003115.1, NM 003115.2, NM 003115.3, NM 003115.4, BC009377.

BC009377.2

UniProt ID: Q16222

Summary: Converts UTP and GlcNAc-1-P into UDP-GlcNAc, and UTP and GalNAc-1-P into UDP-GalNAc.

Isoform AGX1 has 2 to 3 times higher activity towards GalNAc-1-P, while isoform AGX2 has 8

times more activity towards GlcNAc-1-P.[UniProtKB/Swiss-Prot Function]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our custom shRNA service.




Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:




GFP signal was observed under microscope at 48 hours after transduction of TL308538A virus into HEK293 cells. TL308538A virus was prepared using lenti-shRNA TL308538A and [TR30037] packaging kit.






GFP signal was observed under microscope at 48 hours after transduction of TL308538B virus into HEK293 cells. TL308538B virus was prepared using lenti-shRNA TL308538B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL308538C] virus into HEK293 cells. [TL308538C] virus was prepared using lenti-shRNA [TL308538C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL308538D] virus into HEK293 cells. [TL308538D] virus was prepared using lenti-shRNA [TL308538D] and [TR30037] packaging kit.