

Product datasheet for **TL308689**

p73 (TP73) Human shRNA Plasmid Kit (Locus ID 7161)

Product data:

Product Type:	shRNA Plasmids
Product Name:	p73 (TP73) Human shRNA Plasmid Kit (Locus ID 7161)
Locus ID:	7161
Synonyms:	P73
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TP73 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 7161). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001126240 , NM_001126241 , NM_001126242 , NM_001204184 , NM_001204185 , NM_001204186 , NM_001204187 , NM_001204188 , NM_001204189 , NM_001204190 , NM_001204191 , NM_001204192 , NM_005427 , NM_005427.1 , NM_005427.2 , NM_005427.3 , NM_001126242.1 , NM_001126242.2 , NM_001126241.1 , NM_001126241.2 , NM_001126240.1 , NM_001126240.2 , NM_001204189.1 , NM_001204186.1 , NM_001204185.1 , NM_001204191.1 , NM_001204184.1 , NM_001204190.1 , NM_001204188.1 , NM_001204187.1 , NM_001204192.1 , BC117251 , BC117253 , NM_001204190.2 , NM_001204186.2 , NM_001204185.2 , NM_001204192.2 , NM_001204184.2 , NM_001126240.3 , NM_001126242.3 , NM_001204191.2 , NM_001204189.2 , NM_001126241.3
UniProt ID:	O15350
Summary:	This gene encodes a member of the p53 family of transcription factors involved in cellular responses to stress and development. It maps to a region on chromosome 1p36 that is frequently deleted in neuroblastoma and other tumors, and thought to contain multiple tumor suppressor genes. The demonstration that this gene is monoallelically expressed (likely from the maternal allele), supports the notion that it is a candidate gene for neuroblastoma. Many transcript variants resulting from alternative splicing and/or use of alternate promoters have been found for this gene, but the biological validity and the full-length nature of some variants have not been determined. [provided by RefSeq, Feb 2011]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).