

## Product datasheet for **TL308782**

### TMC6 Human shRNA Plasmid Kit (Locus ID 11322)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	TMC6 Human shRNA Plasmid Kit (Locus ID 11322)
Locus ID:	11322
Synonyms:	EV1; EVER1; EVIN1; LAK-4P; Inc; TNRC6C-AS1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TMC6 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11322). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001127198</a> , <a href="#">NM_001321185</a> , <a href="#">NM_007267</a> , <a href="#">NM_007267.1</a> , <a href="#">NM_007267.2</a> , <a href="#">NM_007267.3</a> , <a href="#">NM_007267.4</a> , <a href="#">NM_007267.5</a> , <a href="#">NM_007267.6</a> , <a href="#">NM_007267.7</a> , <a href="#">NM_001127198.1</a> , <a href="#">NM_001127198.2</a> , <a href="#">BC023597</a> , <a href="#">BC023597.2</a> , <a href="#">BC007766</a> , <a href="#">BC018346</a> , <a href="#">BC035648</a> , <a href="#">NM_001127198.5</a>
UniProt ID:	<a href="#">Q7Z403</a>
Summary:	Epidermodysplasia verruciformis (EV) is an autosomal recessive dermatosis characterized by abnormal susceptibility to human papillomaviruses (HPVs) and a high rate of progression to squamous cell carcinoma on sun-exposed skin. EV is caused by mutations in either of two adjacent genes located on chromosome 17q25.3. Both of these genes encode integral membrane proteins that localize to the endoplasmic reticulum and are predicted to form transmembrane channels. This gene encodes a transmembrane channel-like protein with 10 transmembrane domains and 2 leucine zipper motifs. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).